



Basket tool for changing and holding baskets

BSKTOOL & BSKTOOL-DK

The *NEW* QLA Basket Tools BSKTOOL & BSKTOOL-DK are for picking up and moving Dissolution Baskets without touching the mesh or rim and can be used with standard & Distek Knurled rim baskets. This makes it easy to fit baskets onto the shaft without risking damage to the delicate wire mesh body, and extending the useful life of the basket.



BSKTOOL
For
Standard Baskets



BSKTOOL-DK
For
Knurled Rim Baskets

Attaching and detaching baskets is simple with either basket tools.

The Basket Tools can be used to insert baskets onto the QLA 12 position basket holder using minimal force. Both QLA Basket Tools are also dishwasher safe.

