


## *Capsule Sinker Specifications*

### *Helix Coil Type*



#### **CAPWHT-02**

|                      |              |
|----------------------|--------------|
| Material             | Music Wire   |
| Finish               | PTFE Coating |
| Inside Length        | .840 in      |
| Inside Diameter      | .370 in      |
| Wire Diameter        | .050 in      |
| Pitch                | .156 in      |
| End Coil Inside Dia. | .175 in      |
| Maximum Width        | .64 in       |
| No. of Coils         | 8.5          |
| Weight               | 2.05 gm      |
| Maximum Capsule Size | 0            |

#### **Ordering Information:**

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at [www.qla-llc.com](http://www.qla-llc.com) to see all our product line.