


Capsule Sinker Specifications Helix Coil Type



CAPWHT-2S

Material	316 L Stainless Steel
Finish	Electropolished
Inside Length	.840 inches
Inside Diameter	.370 inches
Wire Diameter	.043 inches
Pitch	.156 inches
End Coil Inside Diameter	.190 inches
Maximum Width	.64 inches
Number of Coils	8.5
Maximum Capsule Size	0

Ordering Information:

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.