


Capsule Sinker Specifications

Helix Coil Type



CAPWHT-LG

Material	Music Wire
Finish	PTFE Coated
Inside Length	1.100 inches
Inside Diameter	.410 inches
Wire Diameter	.043 inches
Pitch	.320 inches
End Coil Inside Diameter	.175 inches
Maximum Width	.62 inches
Number of Coils	6.5
Maximum Capsule Size	000

Ordering Information:

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at www.qla-llc.com to see all our product line.