


## *Capsule Sinker Specifications*

### *Helix Coil Type*



#### CAPWHT-LGMT

Material	Music Wire
Finish	Orange Powder Coat
Inside Length	1.100 in
Inside Diameter	.410 in
Wire Diameter	.043 in
Pitch	.320 in
End Coil Inside Dia.	.175 in
Maximum Width	.62 in
No. of Coils	6.50
Weight	2.00 gm
Maximum Capsule Size	000

#### **Ordering Information:**

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your Quality Lab Accessories Sales Representative to order others sizes we carry, Or refer to our web site at [www.qla-llc.com](http://www.qla-llc.com) to see all our product line.