


# ***NEW*** USP APP. 4 Flow Cells *For Flow Through Dissolution Tests*



QLA's flow cells are made of the highest quality to meet or exceed USP Specifications. All flow cells include O-Rings and are available in 12mm and 22.6mm inside diameter sizes and in either Clear Acrylic or Ultem plastic.

Ultem is a semi-translucent amber colored plastic that resists heat to 200 deg. C, and is superior to nearly all other amorphous plastics. QLA's flow cells fit all Sotax CE 7smart Flow Through Dissolution Testers.

**Ordering Information:**

QLA Part No.	Size (mm)	Material	OEM	OEM Cross Ref.
APP4FC-12S	12	Acrylic	Sotax	8386
APP4FC-22S	22.6	Acrylic	Sotax	8220
APP4FC-12SU	12	Ultem	QLA	N/A
APP4FC-22SU	22.6	Ultem	QLA	N/A

03/02/2010