

## Quality Lab Accessories 2 Ton Universal Press

-  Includes a Solid Aluminum Base Stand
-  Simple and Safe 1 Person Operation
-  Incremental Screw Pressure System

The New *QLA 2 Ton Universal Press* utilizes a simple and effective screw system to apply up to 2 Tons of force onto the pellet die. *Gradual Screw Twist Force* causes *Slow and Measured Pressure* to be applied, thereby eliminating uncontrolled pressure jumps that may over pressurize the pellet with other types of systems.



The New *QLA 2 Ton Universal Press* is designed to work with most common Intrinsic Dissolution Apparatus Dies, “Woods Apparatus Design” and specifically the Quality Lab Accessories P/N IDASET-01 Intrinsic Dissolution Apparatus. But will also fit other 4” Round & Square Die Plates made by many manufacturers.

The New *QLA 2 Ton Universal Press* sits on top of a solid ½ inch thick Aircraft Aluminum base plate with pockets for dies, screws and the steel twist bar to fit smartly against the Tablet Press Frame. The base and Press are portable and can be moved from one table location to another. The base plate is secured with 4 steel pins that mount ½ inch into any laboratory tabletop surface, and can be removed by hand.

The entire *QLA 2 Ton Universal Press* frame and base plate are completely Powder Coated Black, and will resist most chemicals and abrasive conditions. Stainless Steel components make for easy cleaning and virtually *Zero Maintenance*.

**Call A Quality Lab Accessories Sales Representative Today To Discuss Your Laboratory Press Needs.**

Part Number	Item
IDAPRS-01	Two Ton Universal Press